

Miracle Belt™

(Techniques & Case Studies)

The Miracle Belt has proved to be a valuable addition to Physical Therapy, Occupational Therapy and Speech Therapy treatment approaches in our clinic. The wearing schedule varies from child to child based on their individual needs and tolerance. The actual application on the belt also varies based on the Child's difficulties. This allows for great versatility in its uses and target populations.

I have had children wear the Miracle Belt for a variety of durations that have included: one hour two to three times per day all the way up to wearing it all morning, taking it off for napping and then wearing it all afternoon. The child will let you know, based on their behavior, when it needs to come off. As with any equipment, the ultimate goal is to wean the child away from needing the belt with carryover of the positive results.

CASE STUDIES:

1) A 2 year 9 month old child with Down Syndrome who is very distractible and aggressive with very limited verbal skills but very good non-verbal skills and many known signs. This child was having a difficult time sitting and participating in Speech Therapy with behaviors such as: pushing the table away, turning around in his chair, swatting at the therapist, pulling off his glasses. The Speech Therapist asked if she could try to use a Miracle Belt and upon placing a 2 pound Miracle Belt around the child's waist with the weight anterior, child immediately sat down, gave the therapist his attention and participated in the remainder of the session without difficulty.

2) A 16 month old little boy with developmental delay, mild left sided weakness, moderate sensory and tactile issues with low frustration tolerance and frequent tantrums. This child began walking about 2 weeks ago and is presenting with a VERY wide base of support, a step to pattern with the right leg, very unstable hips and frequent falling. While in therapy this child had several tantrums when his play was redirected. I tried a 1 pound Miracle Belt with the weight anterior and the effect was initial irritation at the application of the belt followed quickly by calming, decreased base of support, increased control, and better balance with less falling. The parent borrowed the Miracle Belt for 2 days and was amazed at how calm, cooperative and steady her child became while wearing the Miracle Belt. The parent also reported that although his base of support is still wide, it was decreased even after removing the Miracle Belt. This child wore the Miracle Belt all day with it being removed for napping and then put back on for the remainder of the day.

3) A 3 year old boy with developmental delays, drug exposure, low tone, fearful of new activities and all imposed movement, significant weakness and resistive to all upper extremity weight bearing activities. I tried the 2 pound Miracle Belt, with the center anterior, during a Physical Therapy session to try to increase focus, increase body awareness and decrease fear. When asked to play in four point (a most hated position) this child complied willingly and was able to weight shift to both the left and right with verbal or tactile cueing without difficulty or complaint. He maintained this position for approximately 2 minutes and completed a puzzle without complaint and in good alignment.

I have seen positive results in many children with a wide variety of needs, but the Miracle Belt is not for everyone. I have had several children that even when the Miracle Belt was re-introduced in later sessions, no benefit was noted. But for many of our children, it has provided an easy, inexpensive alternative that provides a very different input than a weighted or compression vest. Because the Miracle Belt is worn around the waist, it provides a greater sense of grounding through the core. The Miracle Belt can be placed with the center of the weight anterior, posterior or over the left or right hip to help the child become more aware of where he (or she) needs to shift his center of gravity. The ability to easily cue a child into their body has proven to be very powerful.

I would suggest trying the Miracle Belt on, varying the weight distribution and assessing the input you perceived. I would also suggest trying it on others that were able to give good verbal feedback so you have a good representation of how it is perceived to many people.

The Miracle Belt has been a valuable addition to our clinic.

Ouida Wellenberger
PT (Physical Therapist)
THE KIDSPOT